



Spare part list 2019

Update May 2019

Nursing care bed practico ultra low 9.5/80


From June 2019 the new generation of the lifting columns are used (maked in red)
The new and old genertion cannot be interchanged with each other!



Please mention the serial number of the bed on your order

Subject to technical modification

Issued 05/2019



Pos.-Nr.	Stück	Benennung	Art.Nr.	Zg.- Nr.
1	1	BG-LF PRACTICO ULTRANIEDRIG	91987	015919
2	1	BG-ENDSTUECK-HS-KOPF-KOMPL.	91964	014578
3	1	BG-ENDSTUECK-HS-FUSS-KOMPL.	91966	014950
4	2	BG-SEITENGITTER-OBJEKT-100X30-Z20+200KOMPL.	91703	013178


Nursing care bed practico ultra low

9.5/80

No.	Article code	Description
1	91987 A044	Complete lying surface with motors and binding beam (removable aluminium slats)
	91987 A043	Complete lying surface with motors and binding beam (wooden slats)
2	92267 A001	Complete head end panel with lifting motor Type 31, R24034 Steamed beech
	92267 A058	Complete head end panel with lifting motor Type 31, R24034 Steamed beech 2. Generation
	92267 A020	Complete head end panel with lifting motor Type 31, R20128 Sonoma Oak light
	92267 A062	Complete head end panel with lifting motor Type 31, R20128 Sonoma Oak light 2. Generation
	92267 A029	Complete head end panel with lifting motor Type 31, R27017 Maple dark
	92267 A063	Complete head end panel with lifting motor Type 31, R27017 Maple dark 2. Generation
	92267 A016	Complete head end panel with lifting motor Type 31, R30011 Walnut dark
	92267 A064	Complete head end panel with lifting motor Type 31, R30011 Walnut dark 2. Generation
	92267 A018	Complete head end panel with lifting motor Type 31, R42009 Cherry planked
	92267 A065	Complete head end panel with lifting motor Type 31, R42009 Cherry planked 2. Generation
	92267 A002	Complete head end panel with lifting motor Type 32, R24034 Steamed beech
	92267 A060	Complete head end panel with lifting motor Type 32, R24034 Steamed beech 2. Generation
	92267 A021	Complete head end panel with lifting motor Type 32, R20128 Sonoma Oak light
	92267 A066	Complete head end panel with lifting motor Type 32, R20128 Sonoma Oak light 2. Generation
	92267 A023	Complete head end panel with lifting motor Type 32, R27017 Maple dark
	92267 A067	Complete head end panel with lifting motor Type 32, R27017 Maple dark 2. Generation
	92267 A004	Complete head end panel with lifting motor Type 32, R30011 Walnut dark
	92267 A059	Complete head end panel with lifting motor Type 32, R30011 Walnut dark 2. Generation
	92267 A003	Complete head end panel with lifting motor Type 32, R42009 Cherry planked
	92267 A068	Complete head end panel with lifting motor Type 32, R42009 Cherry planked 2. Generation
3	91966A001	Complete foot end panel with lifting motor Type 31, R24034 Steamed beech
	91966 A080	Complete foot end panel with lifting motor Type 31, R24034 Steamed beech 2. Generation

	91966 A038	Complete foot end panel with lifting motor Type 31, R20128 Sonoma Oak light
	91966 A084	Complete foot end panel with lifting motor Type 31, R20128 Sonoma Oak light 2. Generation
	91966 A049	Complete foot end panel with lifting motor Type 31, R27017 Maple dark
	91966 A085	Complete foot end panel with lifting motor Type 31, R27017 Maple dark 2. Generation
	91966 A050	Complete foot end panel with lifting motor Type 31, R30011 Walnut dark
	91966 A086	Complete foot end panel with lifting motor Type 31, R30011 Walnut dark 2. Generation
	91966 A017	Complete foot end panel with lifting motor Type 31, R42009 Cherry planked
	91966 A087	Complete foot end panel with lifting motor Type 31, R42009 Cherry planked 2. Generation
	91966 A006	Complete foot end panel with lifting motor Type 32, R24034 Steamed beech
	91966 A082	Complete foot end panel with lifting motor Type 32, R24034 Steamed beech 2. Generation
	91966 A051	Complete foot end panel with lifting motor Type 32, R20128 Sonoma Oak light
	91966 A088	Complete foot end panel with lifting motor Type 32, R20128 Sonoma Oak light 2. Generation
	91966 A052	Complete foot end panel with lifting motor Type 32, R27017 Maple dark
	91966 A089	Complete foot end panel with lifting motor Type 32, R27017 Maple dark 2. Generation
	91966 A047	Complete foot end panel with lifting motor Type 32, R30011 Walnut dark
	91966 A081	Complete foot end panel with lifting motor Type 32, R30011 Walnut dark 2. Generation
	91966 A025	Complete foot end panel with lifting motor Type 32, R42009 Cherry planked
	91966 A090	Complete foot end panel with lifting motor Type 32, R42009 Cherry planked 2. Generation
		OPTIONAL :
4	91703 A001	Pair of continuous wooden side rails (extendable)for one bed side R24034 Steamed beech
	91703 A058	Pair of continuous wooden side rails (extendable)for one bed side R20128 Sonoma Oak light
	91703 A061	Pair of continuous wooden side rails (extendable)for one bed side R27017 Maple dark
	91703 A050	Pair of continuous wooden side rails (extendable)for one bed side R30011 Walnut dark
	91703 A004	Pair of continuous wooden side rails (extendable)for one bed side R42009 Cherry planked


Liegefläche practico ultraniedrig 9,5/80	91987
--	-------



Pos-Nr	Stück	Benennung	Art-Nr.	Zg-Nr
1	1	BG-LIEGEFLÄCHE PRACTICO ULTRA NIEDRIG UNTERZUG	9195792010	0157140
2	1	BG-PRACTICO ULTRA NIEDRIG UNTERZUG	92181	015906
3	4	DIN912-M8X35	700.00545	003225
4	2	Gelenkschraube-M8X45	700.00548	000955
5	4	U-Scheibe 84X25X3, verzinkt	700.00466	014949
6	6	Stopmutter-M8	700.00446	000731
7	1	MCL-I	800.00635	014490
8	4	Senkschraube M4X12	700.00623	007146
9	4	Blind-Einnetmutter-M4	700.00985	015214
10	4	Rückwandschraube 5X35-Kreuz	700.00926	014848
11	1	Befestigungsschnelle 5MM	890.01479	015948
12	1	Bohrschraube 4,8X16-M-Kreuz	700.00799	000446

Lying surface practico ultra low 90x200 cm

No.	Article code	Description
1	91957 A003	Lying surface incl. Motors (Removable aluminium bars)
	91957 A002	Lying surface incl.motors (Wooden slats)
2	92181 A002	Binding beam
3	700.00545	Cylinder head screw M8x35
4	700.00548	Cylinder head screw M8x45
5	700.00466	Felt washer 8,4x25x3
6	700.00446	Stop nut M8
7	800.00635	Control box Dewert MCL 2
	800.00785	Control box Dewert MCL 3 - starting June 2019
8	700.00623	Countersunk screw M4x12
9	700.00985	Blind rivet M4
10	700.00926	Screw 5,0x35
11	890.01479	Fastening clip 5mm
12	700.00709	Drilling screw filister head 4,8x16




Pos.-Nr.	Stück	Benennung	Art.Nr.	Zg.-Nr.
1	1	BG-ALUMINIUMRAHMEN-PRACTICOALUPLUS-MOD2015	92173	015741
2	1	BG-KOPFHEBETEI.LALUPLUS-KOMPL.	92174	015742
3	1	BG-LASCHE.KB.PRACTICOALUPLUS.LUECKENSCHLUSS LINKS	920.04.389	015771
3	1	BG-LASCHE.KB.PRACTICOALUPLUS.LUECKENSCHLUSS	920.04.389	015770
4	1	BG-FUSSHEBETEI.PRACTICOALUPLUS-KOMPL.	92176	015745
5	1	BG-KB-ALUPLUS-860X790-STAHL-ALU-FEDER-KOMPL.	91953/91954	015746
6	2	WINKEL.KB.HS	610.00703	014589
7	1	BG-MITTELTEI.LOSE-ALUPLUS-KOMP.	92175	015755
8	1	BG-FB-ALUPLUS-500X790-STAHL-ALU-FEDER-KOMP.	91955/91956	015753
9	1	BG-MEGAMAT2-377-160MM	800.00636/800.00676	015997
10	1	MEGAMAT-XS	800.00637/800.00677	014573
11	3	SPLINTBOLZEN-MOTORBEFESTIGUNG-10X50X45	703.00114	012590
12	3	SICHERUNGSLIP-FUER-SPLINTBOLZEN	700.00681	001256
13	1	BG-ROHRKLAPPSTECKER	700.00973	014700
14	4	BG-LEISTENAUFNAHME-ALUSTAHL	91655, 91656, 91657	015748
15	13	BG-ALUROHR-70X14-1FEDER	91655, 91656, 91657	015768
16	4	BG-BUEGEL-RIPOLUX-MIT-SCHALE	92177	015534
17	1	BG-ENSCHUB-ABNEHMBENDST-RE-MOD2014-KOMP.	920.04.165	013624
18	1	BG-ENSCHUB-ABNEHMBENDST-LI-MOD2014-KOMP.	920.04.166	013622
19	1	BG-ENSCHUB-ABNEHMBAR-ENDSTUECK-MOD2014-RE	920.04.164	013164
20	1	BG-ENSCHUB-ABNEHMBAR-ENDSTUECK-MOD2014-LI	920.04.163	013163
21	4	BG-BEFESTIGUNG.ENSCHUB.ALUBETT	910.084.84	015210

Lying surface practico ultra low


No.	Article code	Description
1	92173 A003	Aluminium frame of the lying for alu lying surface
	92173 A002	Aluminium frame of the lying for wooden lying surface
2	92174A001	Head bar actuator
3	920.04389	Connecting link right/left
4	92176A001	Foot bar actuator
5	91954A001	Head bar for lying surface with removable aluminium bars
	91953 A001	Head bar for lying surface with wooden slats
6	610.00703	Angle bracket mattress compensation
7	92175A001	Middle part right/left removable aluminium bars
	92175 A002	Middle part right/left wooden slats
8	91956A001	Foot bar for lying surface with removeable alu bars
	91955 A001	Foot bar for lying surface with wooden slats
9	800.00636	Head part single drive Dewert
10	800.00637	Foot part single drive Dewert
	800.00758	Hand control Dewert
	800.00662	Nurse key hand control
11	703.00114	Bolt with nut 10/45,8/52
12	700.00681	Piston bolt SL 10mm
13	700.00973	Pipe linch pin
14	830.00467	Fixation clip aluminium bars
	890.01965	Rubber mould for wooden slats
	700.00815	Metal screw filister head (4.8x16) for wooden slats
15	91345A001	Removable aluminium bar (lying surface)
	91667 A001	Wooden slats
	91669 A001	Wooden slats with tension regulators
16	92177A001	Mattress bracket
17	920.04165	End panel fixation foot end, right side
18	920.04166	End panel fixation foot end, left side
19	920.04164	End panel fixation head end right side
20	920.04163	End panel fixation head end left side
21	910.08484	Connecting link alu bed

Endstück kopfseitig	91964
---------------------	-------



Pos.-Nr.	Stück	Benennung	Art.Nr.	Zg.-Nr.
1	1	BG-ENDSTUECK-HS-KOPF-MONTIERT	91989	015912
2	1	BG-ENDSTUECK-TYP-32-HS	91961	014577
2	1	BG-ENDSTUECK-TYP31-KOMP	91961	014969
2	1	ENDSTUECK-TYP33-OHNE_SG	80424	014986
3	4	FLACH-KOPFSCHRAUBE-M6X30	70000043	001467


Endstück kopfseitig vormontiert	91989
---------------------------------	-------



Pos.-Nr.	Stück	Benennung	Art.Nr.	Zg.-Nr.
1	1	BG-HS-T14	80000649	014566
2	2	BLIND-ENNETMUTTER-M6	70000051	009202
3	1	HPL-DEKOR-GEBOGEN	80401	014834
4	2	GEWINDEPLATTE_HS_3XMB_1X16	61000712/91008589	014697
5	1	BG-ANBINDUNG_HS_KOPF	91962	014583
6	8	ZYLINDERKOPFSCHRAUBE-M8X20	70000539	014993
7	1	BG-FUSSROHR-HUBSAEULE_KMPL	92182	014845
8	4	DIN12-M8X70	70000556	000965
9	1	HAUBE_HUBSAEULE	83000528	014107
10	1	ZUGENTLASTUNG_ANBINDUNG_HS	61000717 (91967)	014942
11	2	ZYLINDERKOPFSCHRAUBE-M8X25	70000542	000963
12	1	PG-9		014943

Endstück fußseitig


91966



Pos.-Nr.	Stück	Benennung	ArtNr.	Zg- Nr.
1	1	BG-ENDSTUECK-HS-FUSS-MONTERT	91944	015921
2	1	BG-ENDSTUECK-TYP-32-HS	91961	014577
2	1	BG-ENDSTUECK-TYP31-KOMP	91961	014969
2	1	ENDSTUECK-TYP33-OHNE_SG	80424	014986
3	4	FLACHKOPFSCHRAUBE-M6X30	700.00043	001467

Endstück kopfseitig vormontiert

91944



Pos.-Nr.	Stück	Benennung	ArtNr.	Zg- Nr.
1	1	BG-HS_1114	800.00649	014566
2	2	BLIND-ENNETMUTTER-M6	700.00051	009202
3	1	HPL-DEKOR-GEBOGEN	80401	014834
4	2	GEWINDEPLATTE_HS_3XMB_10X6	610.00712/910.08589	014697
5	1	BG-ANBINDUNG_HS_FUSS	91963	014946
6	10	ZYLINDERKOPFSCHRAUBE-M6X20	700.00539	014993
7	1	BG-FUSSROHR-HUBSAEULE_KMPL	92162	014845
8	4	DN972-M6X70	700.00556	000965
9	1	HAUBE_HUBSAEULE	830.00528	014107
10	2	ZYLINDERKOPFSCHRAUBE-M6X25	700.00542	000963

Head end panel


No.	Article no.	Description
1	91989 A001	Lifting coloumn ultra low head end
	92266 A002	Lifting coloumn ultra low head end 2. generation
2	92224 A003	End panel type 31, R24034 Steamed beech
		<i>alternativ:</i>
	92224 A001	End panel type 32, R24034 Steamed beech
	92224 A030	End panel type 31, R20128 Sonoma Oak light
		<i>alternativ:</i>
	92224 A031	End panel type 32, R20128 Sonoma Oak light
	92224 A041	End panel type 31, R27017 Maple dark
		<i>alternativ:</i>
	92224 A033	End panel type 32, R27017 Maple dark
	92224 A038	End panel type 31, R30011 Walnut dark
		<i>alternativ:</i>
	92224 A027	End panel type 32, R30011 Walnut dark
	92224 A023	End panel type 31, R42009 Cherry planked
		<i>alternativ:</i>
	92224 A010	End panel type 32, R42009 Cherry planked
3	700.01018	Cylinder head screw M6x30

Foot end panel

No.	Article no.	Description
1	91944 A001	Lifting coloumn ultra low foot end
	91944 A002	Lifting coloumn ultra low foot end 2.generation
2	91961A001	End panel type 31, R24034 Steamed beech
		<i>alternativ:</i>
	91961A006	End panel type 32, R24034 Steamed beech
	91961 A038	End panel type 31, R20128 Sonoma Oak light
		<i>alternativ:</i>
	91961 A049	End panel type 32, R20128 Sonoma Oak light
	91961 A051	End panel type 31, R27017 Maple dark
		<i>alternativ:</i>
	91961 A052	End panel type 32, R27017 Maple dark
	91961 A054	End panel type 31, R30011 Walnut dark
		<i>alternativ:</i>
	91961 A047	End panel type 32, R30011 Walnut dark
	91961 A055	End panel type 33, R30011 Walnut dark
	91961 A018	End panel type 31, R42009 Cherry planked
		<i>alternativ:</i>
	91961 A029	End panel type 32, R42009 Cherry planked
3	700.01018	Cylinder head screw M6x30

Fußrohr komplett

92182




Pos.-Nr.	Stück	Benennung	ArtNr.	Zg.- Nr.
1	1	BG-FUSSROHR-HUBSÄULE	91943	014844
2	2	ZYLINDERSCHRAUBE-M10X90	700.01001	015946
3	2	U-SCH-EBE-20X2-D=105-DIN125	700.00472	000635
4	2	DLR-5947UAP100P30-1LKICKPEDAL	890.02157	014815
5	2	WANDABWESROLLE ULTRANIEDRIG	830.00550	014799
6	2	SECHSKANTMUTTER-M10	700.00212	002903
7	2	ROHRVERBINDUNG-15X25-19,5	609.00118	014810
8	1	FEDERFUSSROHR-LS	830.00553	014823
9	1	NUT_FUSSROHR-LS	830.00552	014824
10	1	BLENDE_FUSSROHR-LS	830.00551	014821
11	2	RUECKWANDSCHRAUBE-5X35-KREUZ	700.00926	014848

Foot tube practico ultra low

No.	Article no.	Description
1	91943A001	Castor fixation
2	700.01001	Cylinder head screw M10x90
3	700.00472	Felt washer 10,5
4	890.02157	Double castor
5	830.00550	Wall protection castors
6	700.00212	Counter sunk screw M10
7	609.00118	Bush tube joint D=15,5x2,5x19,5
8	830.00553	Foot cover
9	830.00552	Foot cover (nut)
10	830.00551	Cover foot tube
11	700.00926	Screw 5,0x35

Endstück Typ 31

91961




Pos.-Nr.	Stück	Benennung	ArtNr.	Zg- Nr.
1	1	ENDSTUECK-HOLZ-MIT BUEGEL-TYP31	804.11	014.970
2	1	ALU-C-SCHENE-F18MM-RECHTS-TYP-31-644_HS	603.00080	014.975
3	1	ALU-C-SCHENE-F18MM-LINKS-TYP-31-644_HS	603.00081	014.974
4	2	BG-SEITENGITTERAUFNAHME-UNIVERSAL	920.04.160	012698
5	2	07981-ST4.8X32-C-H	700.00026	013205
6	4	SCHRAUBE-M6X20-KREUZ-DIN 7985	700.00177	000842
7	6	HOLZSCHRAUBE-4X20-KREUZ-DIN 7505	700.00373	000847
8	1	STOPFEN-ENDSTUECK-2012-UNTEN-LI	830.00471	01B16
9	1	STOPFEN-OBEN-LI-ENDSTUECK	830.00507	012833
10	1	STOPFEN-ENDSTUECK-2012-UNTEN-RE	830.00472	01B18
11	1	STOPFEN-OBEN-RE-ENDSTUECK	830.00506	012831
12	2	STOPFEN-OBEN-INNEN-ENDSTUECK-90ER	830.00499	013033
13	4	RUNDMUTTER-RAMPA-M6-D=17X14,5	701.00024	000844
14	2	SPANNPLATTENSCHRAUBE-4,0X35	700.00375	004.172
15	2	STOPFEN-ENDSTUECK-INNEN-MOD.2014	830.00512	013360
16	2	HUELSFENDSTUECK-20X10X25MM	704.00006	000435
17	1	ROHRBUEGEL-F-ENDSTUECK-TYP-21-MOD.2014	91701	012916

Wooden end panel type 31

No.	Article No.	Description
1		Not available as a spare part
2	603.00080	C-bar side rails type 31 right side
3	603.00081	C-bar side rails type 31 left side side
4	920.04160	Universal side rail fittings (kit)
5	700.00026	Sheet metal screw filister head 4,8x32
6	700.00177	Fillister head screw M8x20
7	700.00373	Chipboard screw with countersunk head 4,0x20
8	830.00471	Cover plug end panel bottom-left side
9	830.00507	Cover plug end panel top-left side
10	830.00472	Cover plug end panel bottom- right side
11	830.00506	Cover plug end panel top-right side
12	830.00499	Cover plug top-inside
13	700.01024	Cylinder head screw M8 D=17x10
14	700.00375	Chipboard screw with countersunk head 4,0x35
15	830.00512	Cover plug end panel inside
16	704.00006	Bush D=10x20 Di=5
17	91701A001	Bow type 31, coated

Endstück Typ 32

91961



Pos.-Nr.	Stück	Benennung	ArtNr.	Zg.- Nr.
1	1	ENDSTUECK-HOLZ-RADIUS-TYP32	804.02	014606
2	1	ALU-C-SCHENE-F19MM-RECHTS-TYP-32-644-JHS	60300079	014609
3	1	ALU-C-SCHENE-F19MM-LINKS-TYP-32-644-JHS	60300078	014608
4	2	BG-SEITENGITTERAUFNAHME-UNIVERSAL	920.04.160	012698
5	2	D7981-ST4.8X32-C-H	700.00026	013205
6	4	SCHRAUBE-M8X20-KREUZ-DIN 7985	700.00177	000842
7	4	HOLZSCHRAUBE-4X20-KREUZ-DIN 7505	700.00373	000847
8	1	STOPFEN-ENDSTUECK-2012-UNTEN-LI	830.004.71	011316
9	1	STOPFEN-ENDSTUECK-2012-OBEN-LI	830.004.73	011218
10	1	STOPFEN-ENDSTUECK-2012-UNTEN-RE	830.004.72	011318
11	1	STOPFEN-ENDSTUECK-2012-OBEN-RE	830.004.74	011317
12	2	STOPFEN-ENDSTUECK-INNEN-MOD.2014	830.00512	013360
13	4	RUNDMUTTER-RAMPA-M8-D=17X14,5	70100024	000844
14	2	SPANNPLATTENSCHRAUBE-4.0X35	700.00375	004172
15	2	HUELSFENDSTUECK-20X10X25MM	704.00006	000435

Wooden end panel type 32

No.	Article No.	Description
1		Not available as as single spare part
2	603.00079	C-bar side rails type 32 right side
3	603.00078	C-bar side rails type 32 left side side
4	920.04160	Universal side rail fittings (kit)
5	700.00026	Sheet metal screw filister head 4,8x32
6	700.00177	Fillister head screw M8x20
7	700.00373	Chipboard screw with countersunk head 4,0x20
8	830.00471	Cover plug end panel bottom-left side
9	830.00473	Cover plug end panel top-left side
10	830.00472	Cover plug end panel bottom- right side
11	830.00474	Cover plug end panel top-right side
12	830.00512	Cover plug top-inside
13	700.01024	Cylinder head screw M8 D= 17x10
14	700.00375	Chipboard screw with countersunk head 4,0x35
15	704.00006	Bush D=10x20 Di=5